

# Institute of Mathematics, Physics and Chemistry

# **Department of Chemistry**

Water chemistry laboratory

Laboratory exercise

Determination of chloride ion content and specific conductivity

Elaborated by: dr inż. Jan Krupowies mgr inż. Czesław Wiznerowicz dr inż. Agnieszka Kalbarczyk-Jedynak dr inż. Konrad Ćwirko dr Magdalena Ślączka-Wilk

KIEROWNIK akładu Chemii 6 Car Ce inz, Agnieszka Kalbarczyk-Jedynak

Szczecin 2023

# **EXERCISE SHEET**

1	Relation to subjects: ESO	/25, 27 DiRMiUO/25,	27 EOUNiE/25, 27								
	Specialty/Subject	Learning outcomes	Detailed learning outcomes								
		for the subject	for the subject								
	ESO/26 Chemistry	EKP3	SEKP6 – Determination of selected								
	of water, fuels	K_U014, K_U015,	indicators of technical water								
	and lubricants	K_U016.	quality;								
	DiRMiUO/26 Chemistry	EKP3	SEKP6 – Determination of selected								
	of water, fuels	K_U014, K_U015,	indicators of technical water								
	and lubricants	K_U016.	quality;								
	EOUNiE/26 Chemistry	EKP3	SEKP6 – Determination of selected								
	of water, fuels	K_U014, K_U015,	indicators of technical water								
	and lubricants	K_U016.	quality.								
2	Purpose of the exercise:										
	0	1 ·	ermine the content of chloride ions								
			learning the technical requirements								
		ods of their measureme	ent and operational significance;								
3	Prerequisites:										
			and safety regulations in a laboratory								
	-	-	e on the appropriate form, knows -								
			of chloride ions and the conductivity								
			and its operating principle, knows								
			ion of water on the ship, as well								
	as technical requirements a		g of marked parameters;								
4	Description of the laborat										
	• •		vsis, salt meter, conductivity meter,								
~	instructions for operating the	he instruments, water s	samples;								
5	Risk assessment *:										
	-	ution of silver nitrate (	0.028M – very low probability of dark								
	spots on the skin. Final assessment – <b>VERY SMALL THREAT</b>										
	Safety measures required	•									
	a. lab coats,	na producta popor tou	vale								
6	b. health and safety cleanies <b>The course of the exercise</b>	** **	ve1s,								
0			d familiarize with the laboratory kit								
	for the exercise,	anuar (appendix 1) an	a familiarize with the faboratory kit								
		of chloride ions us	ing the Mohr method and measure								
			measurements with the conductivity								
	meter of the tested wate	-	measurements with the conductivity								
7	Exercise report:										
/		accordance with the in	structions contained in the workplace								
	manual (appendix 1),		structions contained in the workplace								
		ined results of measure	ements and determinations, determine								
			f the tested water by comparing								
	the determined paramet	-	, , ,								
	-	-	t or appropriate corrective action,								
			estions included in the set of tasks								
	and questions to be con	-									
L		ipieces of the student,									

8	Archiving of research results:
	Submit a written report on the performed exercise to the teacher.
9	Assessment method and criteria:
	a. EKP1, EKP2 – tasks given for independent solution and development:
	mark 2.0 – has no basic chemical and operational knowledge concerning
	the determined operational parameters of the tested boiler or cooling water, i.e.
	the content of chloride ions, salinity and specific conductivity of the water; mark 3.0 – has basic chemical and operational knowledge concerning
	the determined functional parameters of the tested boiler or cooling water, i.e.
	the content of chloride ions, salinity and specific conductivity of water,
	and the ability to make basic chemical calculations and solve simple tasks
	in the field of these parameters;
	mark $3.5 - 4.0$ – has extended chemical and operational knowledge in the field
	of determined utility parameters of the tested water and the ability to solve
	complex tasks in the field of assessing changes in these parameters;
	mark $4.5 - 5.0$ – has the ability to apply complex chemical and operational
	knowledge to partial evaluation of the quality and operational suitability of the tested water due to the content of chloride ions, salinity and specific
	conductivity of water, and the ability to make diagnostic decisions on this basis,
	as well as corrective and remedial actions.
	b. EKP3 – control works:
	mark $2.0 - does$ not have the ability to analyse and evaluate the results
	of the performed analyses and determinations and to draw conclusions;
	mark $3.0$ – has the ability to analyse the obtained results, interpret the laws
	and phenomena, transform formulas, and interpret charts and tables; mark $2.5 - 4.0$ , has the shility to breaden the analysis of results, apply laws
	mark $3.5 - 4.0$ – has the ability to broaden the analysis of results, apply laws, construct monograms and charts;
	mark $4.5 - 5.0$ – has the ability to comprehensively analyse the obtained results
	of determinations and measurements, make generalizations, detect cause-and-
	effect relationships and make appropriate operational decisions.
10	References:
	1. J. Krupowies, C. Wiznerowicz, A. Kalbarczyk-Jedynak, K. Ćwirko, M. Ślaczka-
	Wilk, Workplace instruction for laboratory exercise: "Oznaczanie zawartości
	jonów chlorkowych oraz przewodnictwa właściwego", 2022 (in Polish).
	2. <u>https://openstax.org/details/books/chemistry-2e (accessed: 10 July 2023)</u> .
	3. <u>https://planm8.io/blog/marine-boiler-water-treatment</u> (accessed: 10 July 2023).
	4. <u>https://marinersgalaxy.com/boiler-water-test-on-ship-name-of-all-tests/</u> (accessed:10 July 2023).
	<ol> <li><u>https://www.imo.org/en/MediaCentre/HotTopics/Pages/Implementing-the-BWM-</u></li> </ol>
	Convention.aspx (accessed: 29 June 2023).
	6. https://www.imo.org/en/GoogleSearch/SearchPosts/Default.aspx?q=water%20tre
	atment%20on%20ships (accessed: 29 June 2023).
	7. <u>https://www.wilhelmsen.com/product-catalogue/products/marine-chemicals/test-</u>
	kits-and-reagents/water-test-kits/test-kit-for-cooltreat-al/ (accessed: 29 June 2023).
1	

## APPENDIX 1 – MANUAL

## **1. SCOPE OF THE EXERCISE**

- getting acquainted with the workplace instructions for the exercise;
- determination of chloride ions in cooling or boiler water using the argentometric method (Mohr's method);
- determination of the specific conductivity of cooling or boiler water with the use of a conductivity meter;
- assessment of the quality of the tested water and suitability for operational purposes and its possible treatment.

#### **2. THEORETICAL INTRODUCTION TO THE EXERCISE**

#### 2.1. Chloride ion content

In terms of chemistry, sea water is a solution of many chemical compounds, mainly various salts, such as: chlorides, bromides, sulphates, carbonates, etc., occurring in various amounts depending on the considered body of water. The most abundant compound in seawater is sodium chloride NaCl (0.9 - 2.6%).

The good solubility of chlorides and their widespread occurrence in the form of natural salt deposits (NaCl, MgCl<sub>2</sub>) make these salts extremely widespread in natural waters.

Chloride ions are the ingredients that induce the taste of the water. The salty taste in its pure form is caused only by sodium chloride. The bitter taste of sea water is caused by the magnesium salts.

Chloride concentration in waters for various purposes is regulated by law. The applicable national regulations limit e.g., the chloride content in drinking water and economic needs at a level not higher than 250 mg/dm<sup>3</sup>. Water containing increased amounts of chloride is corrosive to concrete, cast iron and steel. For example, groundwater shows a destructive effect even at the content of 100 mg  $Cl^{-}/dm^{3}$ , especially in the case of soft waters and waters containing magnesium chloride.

Sea water, due to its high salinity, cannot be used as feed water for marine boilers. These boilers must be supplied with water of appropriate quality, the so-called fresh water, which is a distillate obtained in ship evaporators by evaporating seawater. Modern distillation installations used on ships produce fresh water from sea water with a salinity ranging from a few to several mg  $Cl^{-}/dm^{3}$ . These devices are equipped with salt meters enabling the continuous control of the content of chloride ions in the produced distillate. The control of the salinity of the distillate can also be carried out using the Mohr method.

The increase in chloride concentration in the distillate is evidence of a malfunctioning evaporator. The change in the chloride content in water is a very important indicator of the technical condition of the water system in use. In a ship's engine room, we deal with mutual contact and interaction of freshwater and sea water circuits (e.g., in coolers, condensers, etc.). In the event of failure of the heat exchangers, there is a risk of seawater entering the freshwater circuit, which is a very harmful phenomenon. The simplest method of detecting such a failure is the systematic control of the concentration of chloride ions in fresh water. The constant increase in the content of chloride ions proves the process of seawater entering the fresh water cycle. In this way, the water analysis carried out on the ship become a kind of diagnostic method for devices operating in water systems.

Thorough control of chloride ions in water is additionally dictated by the fact that chemical water treatment preparations used in shipbuilding do not reduce chloride content, but only remove hardness (descale) and prevent corrosion and water foaming. The permissible limit concentration of chloride ions depends on the type of boiler and its operating parameters. For some types of marine boilers, they are given in the auxiliary tables at the end of the manual.

Excessive salinity of boiler water causes not only increase of its electrolytic conductivity, density and tendency to foaming, but also intensification of corrosive properties. For example, excessive chloride content in water creates a risk of internal boiler corrosion due to their hydrolytic decomposition into aggressive hydrochloric acid. For example, magnesium chloride undergoes hydrolytic degradation at elevated temperatures according to the following reaction:

 $MgCl_2 + 2H_2O \rightarrow Mg(OH)_2 + 2HCl$ 

The hydrochloric acid thus formed is then corrosive to the metal components of the boiler according to the reaction:

 $Fe + 2HCl \rightarrow FeCl_2 + \uparrow H_2$ 

Magnesium hydroxide as sparingly soluble, precipitates in the form of a sediment, increasing the amount of silt and limescale. Chlorides can also release hydrochloric acid by reacting with sulfates in the boiler water, according to the reaction:

 $2NaCl + MgSO_4 + 2H_2O \rightarrow Na_2SO_4 + \downarrow Mg(OH)_2 + 2HCl$ 

#### 2.2. Conductivity

The flow of electric current in the material environment consists in the transport of electric charge, the carriers of which may be electrons or ions. Due to the mechanism of current conduction, conductors are divided into: electronic or metallic conductors (metals, semiconductors, superconductors) and ion or electrolytic conductors. Typical ion conductors are solutions of salts, acids and bases, molten or solid salts.


In the case of electrolytic conductors, the current flow is related to the movement of ions. This causes chemical changes and changes in the concentration of conductive electrolytes. As a rule, electrolytes are much worse conductors of electricity than metals. Their specific conductivity is less than 1S/cm, while the best metallic conductors have a conductivity of 104 - 106 S/cm. Generally, the specific resistance of electrolytes decreases with increasing temperature. The conductivity of electrolyte solutions depends on their concentration, temperature, chemical nature of the solute (ability to dissociate into ions) and the type of solvent (ability to form solvates with the formed ions and the value of electric permeability).

The flow of current through the electrolyte solution (electrolytic conductor) is accompanied by chemical changes taking place on the contact surfaces of two types of conductors (discharge and ion release on the electrodes, secondary reactions with the electrode or solvent) and changes in the concentration of the electrolyte in the spaces near the electrodes. The chemical processes taking place in the analysed systems depend on the chemical properties of the metallic electrodes, the components of the solution and the difference in potential between the electrodes.

Conductometry is an electroanalytical method based on the study of the electrical conductivity of a solution between two electrodes. Testing the conductivity of solutions concerns primarily electrolyte solutions, this conductivity is called electrolytic conductivity or electrolytic conductance.

Direct conductometry consists in measuring the specific conductivity of electrolytes, which for electrolyte solutions is in the range  $0.1 - 10^6 \,\mu\text{S/cm}$ . The determination of the concentration of solutions on the basis of the conductivity test can be performed only for simple electrolyte systems for which there is a linear relationship between conductivity and concentration. The conductivity of most electrolyte solutions does not differ much (except H<sup>+</sup> and OH<sup>-</sup>), therefore conductometry is a non-selective method.

Direct conductometry has been used, among others, in in the industrial analysis of solutions, both concentrated and diluted. It is often used to control technological processes, including water purification technology. It serves, among others to characterize distilled water, boiler feed water, to characterize wastewater.



The values of specific conductivity of some solutions and for water are shown in Fig. 1.

Fig. 1. Conductivity of water and some water solutions (source: https://andyjconnelly.wordpress.com/2017/07/14/conductivity-of-a-solution/) (18.10.2021)

Conductivity is a measure of the ability of a solution to conduct an electrolytic current. Pure water has very little conductivity; it grows rapidly after dissolving acids, hydroxides and salts in it (table 1). The increase in specific conductivity, always observed in the area of low concentrations, with increasing concentration is the result of an increase in the number of ions per unit volume of solution. In the range of higher concentrations, the following factors play a dominant role: interionic interactions, which reduce the mobility of ions, and changes in the value of the dissociation constant. The dependence of the conductivity on the concentration provides valuable information about the changes undergone by the electrolyte.

#### Table 1

	Molarity of the electrolyte											
Floatrolyta	0.00	01 M	0.00	02 M	0.00	05 M	0.00	10 M				
Electrolyte	]	Electrolyte concentration – column 1, conductivity – column 2										
	mg/dm <sup>3</sup>	μS/cm	mg/dm <sup>3</sup>	µS/cm	mg/dm <sup>3</sup>	μS/cm	mg/dm <sup>3</sup>	µS/cm				
Formula	1	2	1	2	1	2	1	2				
NaNO <sub>3</sub> (sodium nitrate)	8.50	10.90	17.0	21.8	42.5	54.20	85.0	107.5				
KNO <sub>3</sub> (potassium nitrate)	10.10	13.10	20.2	26.2	50.5	64.95	101.0	129.0				
KCl (potassium chloride)	7.40	13.50	14.9	26.9	37.3	68.00	74.5	137.0				
NaCl (sodium chloride)	5.85	11.30	11.7	22.8	29.2	56.25	58.5	118.0				
KOH (potassium hydroxide)	5.60	25.00	11.2	50.0	28.0	123.00	56.1	244.0				
NaOH (sodium hydroxide)	4.00	21.60	8.0	43.2	20.0	108.00	40.0	216.0				

# Conductivity of electrolytes at 20°C

#### **3. PERFORMING THE EXERCISE**

Fig. 2 shows a general view of the laboratory workplace for testing the content of chloride ions in water.



Fig. 2. Laboratory workplace for determination of chloride ions by Mohr's method

# **3.1.** Determination of chloride ions in boiler or cooling water by Mohr's method

The determination principle is based on the titration of chlorides with silver nitrate against potassium chromate ( $K_2CrO_4$ ) as an indicator in a neutral or slightly alkaline sample (pH 6.5 – 10). Silver ions form sparingly soluble silver chloride AgCl with chloride ions. In the presence of chromates, silver nitrate AgNO<sub>3</sub> first precipitates AgCl (white cheese-like) of silver chloride:

and after complete precipitation of the chloride ions, a reddish-brown silver chromate a  $(Ag_2CrO_4)$  precipitate is formed. The colour change of the titrated water sample indicates that the equivalence point has been reached, i.e., the end of the chloride ion precipitation reaction in the form of AgCl precipitate. Since there are no Cl<sup>-</sup> ions in the solution, the further addition of AgNO<sub>3</sub> reacts with K<sub>2</sub>CrO<sub>4</sub> according to the reaction:

The order of precipitation of AgCl and Ag<sub>2</sub>CrO<sub>4</sub> sediments (selectivity) is determined by the values of the solubility products of these compounds ( $K_{sp AgCl} < K_{sp Ag_2CrO_4}$ ). Since this value for AgCl is lower, it is less soluble and its critical value is exceeded earlier. So, practically, as long as all chloride ions are not precipitated in the form of AgCl, no Ag<sub>2</sub>CrO<sub>4</sub> precipitate is formed – i.e., it is formed and dissolves immediately. This is the selectivity of the precipitation of the above-mentioned sparingly soluble deposits. The determination is disturbed by acids, hydroxides colour above  $30^{\circ}$  on the platinumcobalt scale. Bromides, iodides and cyanides precipitate like chlorides. Orthophosphates interfere above  $25 \text{ mg/dm}^3$ , precipitating silver phosphate. Iron content above  $10 \text{ mg/dm}^3$  masks the end point of the titration. Sulphides, thiosulfates and sulfites are also disturbing.

They can be eliminated as follows: make the sample alkaline against phenolphthalein with sodium hydroxide solution, add 1 cm<sup>3</sup> of  $H_2O_2$ , mix, then neutralize with sulphuric acid.

If the water has a high colour, add 3  $\text{cm}^3$  of the Al(OH)<sub>3</sub> suspension to 100  $\text{cm}^3$  of the sample, mix, allow the sediment to settle, filter, wash and use the filtrate combined with the washing water for testing.

#### Performing the determination

Pipette 25 cm<sup>3</sup> of the water under test into 3 conical flasks. Then add 75 cm<sup>3</sup> of distilled water to each of them. If interfering factors need to be removed, proceed as outlined above.

Perform the determination in the range of pH 6.5 - 10. If the pH of the sample is outside this range, it must be adjusted with sulphuric acid or sodium hydroxide. Then add 1 cm<sup>3</sup> of 10% K<sub>2</sub>CrO<sub>4</sub> solution as an indicator (Fig. 3). Titrate with standard AgNO<sub>3</sub> solution. The intermediate stage of this titration is shown in Fig. 4, where the turbidity of the solution due to the formation of a white AgCl precipitate is clearly visible. Continue the titration until the colour changes from greenish-yellow to red-brown (Figures 5 and 6). In the same way, titrate the control sample containing 75 cm<sup>3</sup> of distilled water (previously used to dilute the test water). Usually, the consumption of AgNO<sub>3</sub> solution per control sample is 0.2 - 0.3 cm<sup>3</sup>.



Fig. 3. Sample prepared for titration



Fig. 5. Final titration effect



Fig. 4. Intermediate stage of titration



Fig. 6. Comparison of the samples after the titration with the initial one

#### Calculation and reporting of results

Calculate the content of chloride ions according to the formula:

$$C_{Cl^{-}} = \frac{(a-b) \cdot M \cdot 35.45 \cdot 1000}{V} \ [mg \ Cl^{-}/dm^{3}]$$

where:

a – average volume of AgNO<sub>3</sub> solution used for titration of the test sample, cm<sup>3</sup>,

b - volume of AgNO<sub>3</sub> solution used to titrate the control sample, cm<sup>3</sup>,

M – the titre of the AgNO<sub>3</sub> solution, M = 0.028 mol/dm<sup>3</sup>,

35.45 - atomic mass of chlorine,

V – volume of the test water sample used for the determination, cm<sup>3</sup>.

The determined content of chloride ions in the tested boiler water should be converted into the equivalent content of sodium chloride, NaCl.

The method of converting the content of chloride ions into the content of NaCl. For example, if the determined content of chloride ions in the boiler water is 120.00 mg/dm<sup>3</sup>, then the conversion of this content into an equivalent amount of NaCl is as follows:

NaCl $\longrightarrow$ Na <sup>+</sup> + 58.44		$H_2O$			
58.44	NaCl	$\longrightarrow$	$Na^+$	+	Cl-
	58.44				35.45
<u>x</u> 1:	X				120,00

$$C_{\text{NaCl}} = \frac{58.44 \cdot 120.00}{35.45} = 197.82 \text{ [mg NaCl/dm^3]}$$

The concentration of 120.00 mg Cl<sup>-</sup>/dm<sup>3</sup> corresponds to the content of 197.83 mg NaCl/dm<sup>3</sup>.

Fig. 7 shows a laboratory workplace for testing water conductivity.



Fig. 7. Laboratory workplace for determination of conductivity

# 3.2. Measurement of conductivity

Measure the specific conductivity with a digital conductivity meter (Fig. 8) as follows:

Calibrate the conductivity meter before running a series of conductivity measurements. Immerse the measuring electrode in the standard solution with a known conductivity value so that the level of the tested liquid is between the marks of min. max. on the electrode (fig. 8). Turn on the instrument, after the conductivity value has stabilized (if necessary), use the calibration screw to correct the conductivity result so that it is consistent with the conductivity content of the standard solution. Turn off the instrument, then rinse the electrode with distilled water.

#### **Measurement:**

- 1) Immerse the measuring electrode in the tested sample so that the level of the tested liquid is between the marks of min. max. on the electrode (fig. 8) and turn on the instrument;
- 2) After the value has stabilized, read the conductivity result on the instrument display.
- 3) Turn off the instrument, then rinse the electrode with distilled water.
- 4) Take the measurement three times.



Fig. 8. Conductivity meter

#### **Elaboration of the results**

With a conductivity in the range of 10  $\mu S/cm$  and below, the results may differ by 10% of the smallest result.

With a conductivity above 10  $\mu$ S/cm, the results may differ by 5% from the lowest result.

The arithmetic mean of at least 3 measurements should be taken as the final result of the determination.

#### **4. DEVELOPMENT OF THE EXERCISES**

- 1. Compare the obtained results of chloride ion content with their acceptable concentrations for the boiler water of the selected type of ship's boiler or for cooling water.
- 2. Compare the obtained conductivity measurement results with the values of the boiler water conductivity of the selected type of ship boiler or cooling water.
- 3. On the basis of the obtained results, assess the quality and suitability of the tested water for operation in selected types of steam boilers.
- 4. The auxiliary tables 2 4 show the technical requirements for utility water for selected types of steam boilers, and appendixes 2 3 show the requirements recommended by Unitor for boiler and cooling water.

#### **5** THE FORM AND CONDITIONS FOR PASSING THE LABORATORY EXERCISE

- 1. Passing the so-called "entry" before starting the exercise.
- 2. Submission of a correct written report on the completed exercise, which should include:
  - short theoretical introduction,
  - operational significance of the measured parameter,
  - processing of the obtained results according to the position manual.
- 3. Final credit for the test at the end of the semester.

## **Additional tasks**

#### I. Examples of tasks with solutions

- 1. How many milligrams of  $Cl^-$  did the salt sample contain when titrated with 35 cm<sup>3</sup> of 0.11 M n AgNO<sub>3</sub> solution?
  - Solution

The product:  $35 \text{ cm}^3 \cdot 0.11 \text{ mmol/cm}^3$  represents the number of millimoles used in the AgNO<sub>3</sub> titration. the number of millimoles of chloride is equivalent to, because:

 $Ag^+ + Cl^- \iff \downarrow AgCl$ 

so in the analyzed sample there were:  $35 \cdot 0.11 = 3.85$  mmol of chlorides. To calculate the result in milligrams, multiply the number of millimoles by the atomic mass of chlorine:

$$3.85 \cdot 35.45 = 136.5 \text{ mg Cl}^{-1}$$

Answer: The sample contained 136.5 mg of Cl<sup>-</sup>ions.

2. Calculate the concentration of magnesium ions  $Mg^{2+}$  and chloride ions  $Cl^{-}$  in a 0.05 M solution of magnesium chloride  $MgCl_2$ , if its dissociation degree  $\alpha = 75\%$  (i.e., 0.75). Solution:

In the ionic state in the solution, the following MgCl<sub>2</sub> concentration exists:

$$\alpha \cdot C_{\text{total}} = 0.75 \cdot 0.05 = 0.0375 \text{ mol/dm}^3$$

The electrolytic dissociation reaction equation has the form:

$$MgCl_2 \stackrel{H_2O}{\longleftarrow} Mg^{2+} + 2Cl^{-}$$

therefore, the concentration of chloride ions in the solution is 2 times higher than the  $Mg^{2+}$  ions, hence the ion concentrations are respectively:

$$[Mg^{2+}] = 0.0375 \text{ mol/dm}^3$$
,  $[Cl^-] = 2 \cdot 0.0375 = 0.0750 \text{ mol/dm}^3$ 

Answer: The concentration of Mg<sup>2+</sup> is 0.0375 mol/dm<sup>3</sup>, and Cl<sup>-</sup> ions is 0.0750 mol/dm<sup>3</sup>.

#### II. Tasks and questions to be completed by the student

#### Tasks

- 1. How many milligrams of chloride are in a 250-milligram sample of salt if 28.5 cm<sup>3</sup> of a 0.0985 M AgNO<sub>3</sub> solution was used to titrate the chloride after dissolving it in water. Answer: 99.5 mg of chlorides.
- 2. 30.15 cm<sup>3</sup> of 0.1112 M AgNO<sub>3</sub> solution were used to titrate a sample with 0.2230 g of salt dissolved in water. How many percent of NaCl did the analysed salt sample contain? Answer: 87.88% of NaCl.
- What was the molar concentration of AgNO<sub>3</sub> solution if 38.25 cm<sup>3</sup> of AgNO<sub>3</sub> solution were used for 250.5 mg of salt sample containing 89.5% KCl? Answer: 0.07857 M of AgNO<sub>3</sub>.
- 4. How many cm<sup>3</sup> of 0.1515 M AgNO<sub>3</sub> solution will be used for titration of chlorides in a solution containing a weight of 0.5340 g of dissolved SrCl<sub>2</sub> Answer: 44.46 cm<sup>3</sup> of 0.1515 M AgNO<sub>3</sub> solution.
- 5. 1.17g of impure sodium chloride NaCl was dissolved in water and titrated with 0.5M AgNO<sub>3</sub> silver nitrate solution. After adding 36 cm<sup>3</sup> of AgNO<sub>3</sub> solution, the precipitation ceased. Calculate the percentage of NaCl in the contaminated salt. Answer: 90% of NaCl.
- 6. The content of sodium chloride NaCl in water is 8.7 g/dm<sup>3</sup>. Calculate the percentage concentration of this solution and the concentration of Cl<sup>-</sup> chloride ions in mg/dm<sup>3</sup>. The density of the solution is d = 1.1 g/cm<sup>3</sup>.

Answer: the percentage concentration of the solution is 0.79%, and the content of  $Cl^{-}$  ions is 5.28 g/dm<sup>3</sup>.

#### Questions

- 1. Report the principle of determination of chloride ions in water by Mohr's method. Write the appropriate equations of chemical reactions in molecular and ionic form.
- 2. What is the process of desalination and demineralization of water?
- 3. How is the distillate (condensate) obtained on sea vessels?
- 4. What is the purpose and chemistry of determining the chloride ion content in boiler water?
- 5. What diagnostic inference can be made on the basis of systematic control of chloride ion content in ship water systems?
- 6. What is the selectivity of the precipitation: AgCl and  $Ag_2CrO_4$  in the determination of chloride ions by the Mohr method?
- 7. What is the corrosive effect of MgCl<sub>2</sub> in boiler water? Justify your answer with appropriate chemical reactions.
- 8. How is the salt content in water measured with a salt meter?

# Auxiliary tables

Type of water Physicochemical quantities Condensate Distillate Feed water Boiler water ≤ 1200 Content of Cl<sup>-</sup> in ppm 12 + 24 $\leq$ Х General hardness in °dH  $\leq 0.084$  $\leq 0.84$ < 0.56 Х Alkalinity p in ppm 150 + 200Х Х Х Content of PO<sub>4</sub><sup>3–</sup> in ppm 2 + 5Х Х х Exponent of the concentration of oxonium  $6.5 + 9.5^{xx}$ Х Х Х ions, pH Oil content in ppm 0 < 3 traces Х Total salt content in ppm  $\leq 3000$ Х х Х

#### Physicochemical values of water for VL 512/10-01 boilers recommended by the manufacturer

<sup>x</sup> Not included in the standard.

<sup>xx</sup> At a temperature of about 20°C.

Table 3

Acceptable ranges of physicochemical values for "Unitor" boilers operating at various pressures
---


Steam	Alkalin	ity		Ions				
pressure in the boiler	р	т	Cl⁻	PO4 <sup>3-</sup>	${\rm SO_4}^{2-}$	$N_2H_4 \cdot 10^{-1}$	pН	
MPa	ppm	ppm	ppm	ppm	ppm	ppm		
0 - 1.75	100 - 150	2p	200	20 - 50	20 - 50	1 - 2	8.3 – 9	
1.75 - 3.1	100 - 150	2p	200	20 - 50	20 - 50	1 - 2	8.3 – 9	
3.1 - 4.2	100 - 150	2p	200	20 - 50	20 - 50	1 - 2	8.3 – 9	
4.2 - 6,0	100 - 120	2p	50	15 - 30	Х	1 - 1.5	8.3 – 9	
6.0 - 8.0	_	_	30	10 - 25	Х	0.5 - 1	8.3 – 9	

<sup>x</sup> Not included in the standard.

Physicochemical qu	uantities	Pressure in MPa									
			1.4		2.4		4.0				
Feed water		average	acceptable	average	acceptable.	average	acceptable				
General hardness °dH		0.02	0.05	0.015	0.02	0.01	0.015				
Content of O <sub>2</sub>	ppm	0.03	0.05	0.02	0.05	0.02	0.03				
Content of Fe	ppm	0.30	_	0.10	0.2	0.05	0.10				
Content of Cu	ppm	_	—	—	_	—	0.01				
pH value at approx.	20°C	8.50	9.50	8.5	9.5	8.5	9.5				
Oil content	ppm		3.00		2		1				
Content of CO <sub>2</sub>	ppm	_	_	_	25	_	25				
Oxidizability of KMnO4	ppm	_	_	_	_	_	20				
Content of NO <sub>2</sub> <sup>-</sup>	ppm	_	_	_	_	_	0.02				
Content of SiO <sub>2</sub>	ppm	_	t	he size is	determined	according	a				
Conductivity	$\mu Scm^{-1}$	—	to	o the man	ufacturer's in	nstructior	18				
Boiler water											
Alkalinity p	ppm	5-15	2 - 20	2 - 8	2 - 10	1-5	2 –7				
Content of SiO <sub>2</sub>	ppm	_	60	_	40	_	35				
Content of P <sub>2</sub> O <sub>5</sub>	ppm	_	_	10	20	5 - 10	10-20				
Conductivity	$\mu Scm^{-1}$	7000	9000	4000	6000	2000	3000				

Water quality for water-tube boilers operating under pressure up to 4 MPa according to P. Orłowski "Steam boilers in industrial power engineering"

#### Unitor boiler water analytical control report form



EDP no 109256

#### Unitor cooling water analytical control report form

UNITOR chemicals c/o Kjemi Service AS P.O. Box 49 3140 Borgheim NORWAY

Year \_\_\_\_

#### COOLING WATERTREATMENT PROGRAMME SHIPBOARD LOG PAD

SPECTRAPAK 309

Ship \_\_\_\_

Flag	Owner

MAIN ENGINE MANUFACTURER

Make up: Shore Distilled Mixed

PRODUCT Dieselguard NB

\_\_\_\_ Month \_ E M A M J J A S C N D

		JACKETS				CAP	1	TONS			PISTONS CAP TONS								
50 PPM MAX		DATE									DATE								
		≥2700	0	0	0	0	0	0	0	0	≥2700	0	0	0	0	0	0	0	- 0
		2520	0	0	0	0	0	0	0	0	2520	0	0	0	0	0	0	0	0
≤		2340	0	0	0	0		0	0	0	2340	0		0	0	0	0	0	
È	02	2160					-				2160		0						0
Ξ	n N		-	0	0	0	0	0	0	0		0	0	0	0	0	0	0	0
2	ppr	1980	0	0	0	0	0	-	0	0	1980	0	0	0	0	-	0	0	0
CHLORIDE LEVEL 50	as	1800	0	0	-	0	0	0	0	0	1800	0	0	0	-	0	0	0	0
	NITRITE as ppm NO <sub>2</sub>	1620	0	0	0	0	0	0	0	0	1620	0	0	0	0	0	0	0	0
	TR	1440	0	0	0	0	0	0	0	0	1440	0	0	0	0	0	0	0	0
1	Z	1260	0	0	0	0	0	0	0	0	1260	0	0	0	0	0	0	0	0
5		1080	0	0	0	0	0	0	0	0	1080	0	0	0	0	0	0	0	0
ŝ		900	0	0	0	0	0	0	0	0	900	0	0	0	0	0	0	0	-
í		≤720	0	0	0	0	0	0	0	0	≤720	0	0	0	0	0	0	0	0
NORMAL CH	CI	≥100	0	0	0	0	0	0	0	0	≥100	0	0	0	0	0	0	0	0
Į	md	80	0	0	0	0	0	0	0	0	80	0	0	0	0	0	0	0	0
	Ш	60	0	0	0	0	0	0	0	0	60	0	0	0	0	0	0	0	0
2	CHLORIDE ppm	40	0	0	0	0	0	0	0	0	40	0	0	0	0	0	0	0	0
-	LO	20	0	0	0	0	0	0	0	0	20	0	0	0	0	0	0	0	0
	S	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
	pН																		
		lguard NB kg																	
		NB liquid Itr											-						
	Make	Make up Itr																	
		FUEL '	VALVES	5	C	CAP	1	TONS			AUX. GENERATORS CAP TONS COMBINED WITH MAIN ENG.								
		DATE									DATE								
		≥2700	0	0	0	0	0	0	0	0	≥2700	0	0	0	0	0	0	0	0
		2520	0	0	0	0	0	0	0	0	2520	0	0	0	0	0	0	0	0
	~	2340	0	0	0	0	0	0	0	0	2340	0	0	0	0	0	0	0	0
	<b>N</b>	2160	0	0	0	0	0	0	0	0	2160	0	0	0	0	0	0	0	0
5	m	1980	0	0	0	0	0	0	0	0	1980	0	0	0	0	0	0	0	0
5	dd s	1800	0	0	0	0	0	0	0	-	1800	0	0	0	0	0	0	0	0
3	NITRITE as ppm NO <sub>2</sub>	1620	0	0	0	-	0	0	-	0	1620	0	0	0	-	0	0	0	0
V	SITE	1440	0	0	0	0	0	0	0	0	1440	0	0	0	0	-	0	0	-
LEVEL 1200-2400 FTM	E	1260	0	0	0	0	0	0	0	0	1260	0	0	0	0	0	0	0	0
	-	1080	0	0	0	0	0	0	0	0	1080	0	0	0	0	0	0	0	6
		900	0	0	0	0	0	0	0	0	900	0	0	0	0	0	0	0	0
		≤720	0	0	0	0	0	0	0	0	≤720	0	0	0	0	0	0	0	0
	75	≥100									≥100								
	m Cl		0	-	0	0	0	0	0	0		0	0	0	0	0	0	0	0
	mdd	80	0	0	0	0	0	0	0	0	80	0	0	0	0	0	0	0	0
	ш	60	0	0	0	0	0	0	0	0	60	0	0	0	0	0	0	0	0
	JIRIC	40	0	0	0,	0	0	0	0	-	40	0	0	0	0	0	0	0	0
1	CHLORID	20	0	0	0	0	0	0	0	0	20	0	0	0	0	0	0	0	0
	Ċ	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
	pН																		
	Diese	lguard NB kg																	
	Rocor	r NB liquid Itr				1													
	Make	up Itr																	
	Comn	nents																	