



**Maritime University of Szczecin**  
**Faculty of Marine Engineering**  
**Department of Physics and Chemistry**



**Physics Laboratory**

## **Laboratory Description**

**Determination of solids density using Nicholson's aerometer**

English version written by  
dr Agata Kowalska

Szczecin 2017

### Objectives:

- Familiarization with mechanical stability conditions – balance of forces and moment of force (torque).
- Determination of density of solids using Archimedes' principle

### Questions and problems to solve:

- Density, weight, specific gravity and mass – the relationship between these quantities.
- Conditions of mechanical stability
- Archimedes' principle
- Determination of solids density using Nicholson's aerometer.

### Short description:

Nicholson's aerometer (see figure 1) is made of a metal float  $P$ ; two pans  $A$  and  $B$  and connecting rod. The metal float and the pan  $A$  are immersed in distilled water. The system is in a state of equilibrium when the sum of forces acting on it is equal to zero.




Fig. 1. Schematic view of Nicholson's aerometer

The pan  $B$  has to be loaded with weights of mass  $m_1$  so that the position of the marker  $K$  is at the level of water surface. The total weight of the aerometer  $Q$  and weight  $Q_1 = m_1g$  is balanced by buoyant force  $W$  acting on submerged parts of the aerometer:

$$Q + Q_1 = W \quad (1)$$

We remove the weights from the pan  $B$  and put on a body of unknown weight  $Q_c$  and a volume  $V_c$ . We load the pan  $B$  with weights of a mass  $m_2$  so that the position of the marker  $K$  was again at the level of water surface. In this case, buoyancy force is balanced by the weights:  $Q$  (of the aerometer),  $Q_2 = m_2g$  (of weights) and  $Q_c$  (of the examined body).

$$Q + Q_2 + Q_c = W \quad (2)$$

Then we remove the body from the pan  $B$  and attach it to the submerged part of the aerometer  $A$ . We have to do it carefully in order to avoid bubbles gathering on the examined body. Now we load the pan  $B$  with mass  $m_3$ . The connecting rod dives again to the level  $K$ . Now, the total

weight of all the elements of the system is balanced by the buoyant force  $W$  acting on submerged parts of aerometer and by the buoyant force  $W_c$  acting on the submerged body:

$$Q + Q_3 + Q_c = W + W_c \quad (3)$$

Using the equations (1)-(3) we can obtain the relation for weight of the examined body and for the buoyant force acting on this body:

$$Q_c = Q_1 - Q_2,$$

$$W_c = Q_3 - Q_2 \quad (4)$$

The weight of the body and the buoyant force acting on it can be determined by using formulas:  $Q_c = m_c g = \rho_c V_c g$  and  $W_c = \rho_w V_c g$ , where  $\rho_c$  and  $\rho_w$  are the examined body density and the density of distilled water, respectively:

$$d = \frac{\rho_c}{\rho_w} = \frac{Q_c}{W_c} = \frac{Q_1 - Q_2}{Q_3 - Q_2} = \frac{m_1 - m_2}{m_3 - m_2} \quad (5)$$

Knowing the array density of distilled water  $\rho_w$  at a given temperature (corresponding to the temperature at the time point when the measurements were performed), we can calculate the body density  $\rho_c$ :

$$\rho_c = d\rho_w \quad (6)$$

Density of a homogenous body can be also determined from the definition of density, after measuring its mass  $m_c$  and volume  $V_c$ :

$$\rho_c = \frac{m_c}{V_c} \quad (7)$$

#### **Literature:**

1. Resnick R., Halliday D., Walker J., *Fundamentals of Physics*, John Wiley & Sons, INC (available editions).